

## Catalog Number

47452-05-010

**Adapter**, for use with low pressure. Barbed connection for 1/16 inch tubing to 1/8 inch tubing. Polypropylene. ARE-Applied Research brand. 10/PK.



### Additional Info

This adapter connects 1/16th inch (1.6 mm) ID tubing to 1/8th inch (3.2 mm) ID tubing when slid over the barbed fittings.

- Maximum Temperature: 100°C
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for bulk pricing.

Click [HERE](#) for a tubing union to connect 1/16th inch (1.6 mm) ID tubing to 1/16th inch (1.6 mm) ID tubing,

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$42.26**

Printed On 2024/05/20