

Catalog Number

47452-10

Adapter, for use with low pressure. Female luer lock to barbed connection for 1/16 inch tubing. Polypropylene. ARE-Applied Research brand. Bulk Price.



Additional Info

This adapter has a female Luer Lock port on one side and a hose barb connector on the other side for 1/16th inch (1.6 mm) ID tubing. *Bulk packaging.*

This adapter is excellent when connecting 1/16th inch tubing to a syringe.

- Maximum Temperature: 100°C.
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for an economical pack of 10.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Printed On 2024/05/20