

## Catalog Number

47452-15

**Adapter**, for use with low pressure. Female luer lock to female luer lock. Polypropylene. ARE-Applied Research brand. Bulk Price.



### Additional Info

This adapter has a female Luer Lock connection port on both sides.

*This adapter is excellent when using a syringe and you want to another syringe or a port with a male Luer Lock or Luer slip connector. Bulk packaging.*

- Maximum Temperature: 100°C.
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for economical pack of 10.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

Printed On 2024/05/20