

Catalog Number

47452-15-010

Adapter, for use with low pressure. Female luer lock to female luer lock.
Polypropylene. ARE-Applied Research brand. 10/PK.



Additional Info

This adapter has a female Luer Lock connection port on both sides.

This adapter is excellent when using a syringe and you want to another syringe or a port with a male Luer Lock or Luer slip connector.

- Maximum Temperature: 100°C.
- Maximum Pressure: 1,000 psig.
- Click [HERE](#) for bulk pricing.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$42.26

Printed On 2024/05/19