

## Catalog Number

49103-FL

**Performance Qualification Solutions, FL.** Set of Solutions for Qualification of HPLC with Fluorescence Detectors. Chemical Solutions Brand. 1 EA.



### Additional Info

This set of NIST Traceable Solutions are Supplemental for use with the HSQ™ or PQ™ Kits for Performance Qualification of HPLC instruments to Qualify Fluorescence (FL) Detectors. Certificate of Analysis is included with each set.

Set Includes:

- 6 x 30ml Linearity Solutions
- 1 x 30ml Diluent Blank Solution

*Solutions will be supplied with at least six months until the expiration date. They are designed to be used promptly and often. Expiration date is two years after the date of manufacture.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,051.95**

Printed On 2024/05/09