


Catalog Number

49307-502R

Tubing, PEEK, Coiled, 0.005 (0.13 mm) ID, 1/16th inch (1.6 mm) OD.
1225 mm length, HPLC grade, solid red, with tubing connectors. ARE-Applied Research brand. 1 EA.



Additional Info

These Flexible, Coiled PEEK Tubes will expand and contract which will save space in places such as behind an instrument and provide a Pulse Dampener Effect for Parabolic Flow Fluids which will result in a quieter base line and also reduce band broadening due to excess tubing lengths.

This coil is 35 mm wide when resting and 200 mm when fully extended, with a total tubing length of 1225 mm, end to end.

Type 4 FlexiTube, color coded: red for 0.005 inch (0.13 mm) ID

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$59.70

Printed On 2024/05/18