


Slurry Reservoir R1

49330-SR3

Part No.

Slurry Reservoir Kit, All-PEEK, for packing 2.1 mm ID columns that are less or equal to 100 mm in length. Reservoir is 150 mm long with a 4.6 mm ID, complete with all parts and accessories needed to pack



Additional Info

This kit is for packing All-PEEK columns 2.1 mm ID with up to and including 100 mm in length.

The Reservoir Body included in this kit is 150 mm long with 4.6 mm ID and has an internal volume of approximately 2.5 ml. This reservoir is equipped with 1/2-20 external threads on both sides. The adapter has an internal 1/4-28 screw thread on one side and internal 1/2-20 screw threads on the other side. The adapter connects the reservoir to the column being packed.

Kit includes everything you need to pack a bioinert All-PEEK column:

- One Slurry Reservoir - the "body" which is Black PEEK.
- One [End Fitting](#) to connect tubing from a pump to the reservoir.
- Two PEEK sealing rings for the bottom and top connection.
- One [Packing Adapter](#) to connect the reservoir to the column.
- One [ETFE sealing ring](#) for the top of the column connection.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$722.25

Printed On 2024/05/20