


## Catalog Number

49331-SR2

**Slurry Reservoir Kit, All-PEEK, for packing 4.6 mm ID columns that are longer than 100 mm in length.** Reservoir is 300 mm long with a 7.5 mm ID, complete with all parts and accessories needed to pack bioinert,



## Additional Info

This kit is for packing All-PEEK columns 4.6 mm ID that are up to 300 mm in length.

The Reservoir Body included in this kit is 300 mm long with 7.5 mm ID and has an internal volume of approximately 13.2 ml. This reservoir is equipped with 5/8-18 external threads on both sides. The adapter has an internal 1/2-20 screw thread on one side and internal 5/8-18 screw threads on the other side. The adapter connects the reservoir to the column being packed.

Kit includes everything you need to pack a bioinert All-PEEK column:

- One Slurry Reservoir - the "body" which is black PEEK.
- One [End Fitting](#) to connect tubing from a pump to the reservoir.
- Two [PEEK sealing ring / frit spacer](#) for the bottom and top connection.
- One [Packing Adapter](#) to connect the reservoir to the column.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$889.44**

Printed On 2024/05/18