


## Catalog Number

49333-SR3

**Slurry Reservoir Kit, All-PEEK, for packing 7.5 mm ID columns that are greater than 100 mm in length.** Reservoir is 600 mm long with a 7.5 mm ID, complete with all parts and accessories needed to pack



## Additional Info

This kit is for packing All-PEEK columns 7.5 mm ID and greater than 100 mm in length.

The Reservoir Body included in this kit is 600 mm long with a 7.5 mm ID and has an internal volume of approximately 26.5 ml. This reservoir is equipped with a 5/8-18 screw thread on both sides. The adapter has an internal 5/8-18 screw thread on one side and internal 5/8-18 screw threads on the other side. The adapter connects the reservoir to the column being packed.

Kit includes everything you need to pack a bioinert All-PEEK column:

- One Slurry Reservoir - the "body" which is black PEEK.
- One [End Fitting](#) to connect tubing from a pump to the reservoir.
- Three [PEEK sealing ring/ frit spacer](#) for the bottom and top connection.
- One [Packing Adapter](#) to connect the reservoir to the column.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$942.94**

Printed On 2024/05/18