

Catalog Number

49420-COIL-4M

Tubing, Stainless Steel, roll, 0.020 inch (0.508 mm) ID, 1/16 inch (1.59 mm) OD, HPLC Grade, 4 meters roll. ARE-Applied Research brand. 1 EA.



Additional Info

This Stainless Steel tubing has many uses in the laboratory but is designed for HPLC. Welded and annealed, they are tested according to EN 10204/3.1B. These Coils are more economical than pre-cut tubing.

- Internal and Outer Dimension Tolerances: +/- 0.05mm.
- Wall Thickness: 0.541 mm

Click [HERE](#) for Tips on Cutting SS Tubing with a hand tool.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$82.80

Printed On 2024/05/19