

## Catalog Number

49590-F116

**Ferrules**, stainless steel, single cone, for 1/16 inch tubing. ARE-Applied Research brand. 10/PK.



### Additional Info

These single cone ferrules are made to be used with our 316L stainless steel tubing connector male nuts when connecting 1/16th inch SS tubing to columns, fittings or detectors in high pressures zones of HPLC systems.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$37.56**

Printed On 2024/05/18