

## Catalog Number

58701-S-K

**Syringe, Chrom. Disposable, 1ml, Luer Slip, Non Sterile. MicroSolv Brand. (10 bags of 100). 1000/PK.**



### Additional Info

These Luer Slip Chrom Syringes are often chosen for chromatography, with syringe filters or a great choice for any situation needing an inert, non-reactive syringe.

Latex free, and does not contain rubber, silicone oil, styrene or DEHP. Manufactured only from laboratory grade polypropylene and polyethylene. These unique plastic syringes are more chemically resistant than rubber-tipped syringes and have a positive safety stop to help prevent accidental spills.

Click [HERE](#) for larger case quantity.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$607.66**

Printed On 2024/05/19