


Catalog Number

58802-30

Pouring Ring, Replacement, Blue, for HPLC Mobile Phase Reservoir Bottles. GL45mm Threads. ARE-Applied Research brand. 1 EA.



Additional Info

Using this Polypropylene Pouring or "Drip Ring" can make decanting liquids from any bottle with a GL45 Thread easier and prevent that last drop from getting away from you.

Friction fit makes this ring stay in the bottle during pouring and capping.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1.93

Printed On 2024/05/18