

Catalog Number

58807-050

Bottle, HPLC Mobile Phase Storage, Glass, Clear, 50ml, Screw Top, GL32 Threads. MicroSolv Brand. 1 EA.



Additional Info

This Clear, Borosilicate Glass Bottle has a 50ml "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Storage Bottle. The small opening and Enamel Graduations and Marking Patch makes this a convenient way to store concentrated Mobile Phase or its Ingredients.

A Blue, Solid Screw Cap without a liner and a polypropylene Drip Ring is included. This Bottle is *Not "Safety Coated"*.

The Temperature range of these glass bottles is -70° to 500°C. The temperature range of the polypropylene cap is 140°C.

Click [HERE](#) for case quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$24.11

Printed On 2024/05/20