

Catalog Number

58807-2L

Bottle, Mobile Phase Reservoir, Glass, Clear, 2L, Screw Top, GL45 Threads. MicroSolv Brand. 1 EA.



Additional Info

One of our Most Popular Items, this Clear, Borosilicate Glass Bottle has a 2L "Total Fill Capacity" and is perfect as an HPLC Mobile Phase Reservoir. The small opening, one liter Capacity and White Enamel Graduations on the outside of the Bottle help you determine the Solvent use and prevent running dry.

This Bottle can be used with [MicroSolv Filter Degassers](#), other Degassers or disposable [Filter Degassers](#). A Solid Screw Cap without a liner and a polypropylene Drip Ring are included. This Bottle is often used with the [Mobile Phase Delivery System](#) and is *Not "Safety Coated"*.

Click [HERE](#) for case quantity.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$80.91

Printed On 2024/05/19