

## Catalog Number

58820-06

**Cap, Mobile Phase Delivery System**, GL45 threads, blue, PE, for HPLC reservoir bottles. A PTFE insert with 3 holes, PEEK plug and O-ring are included. ARE-Applied Research brand. 1 EA.



### Additional Info

Most labs use this type of bottle cap system (also called solvent delivery caps) to ensure that their mobile phase is not contaminated and thereby helps to comply with their HPLC GLP requirements. Your mobile phase will only come in contact with medical grade PTFE so there will NOT be any contamination from plastic or wax film.

Keeping these polyethylene screw caps with the PTFE center plug on the bottles with the tubing through-holes prevents dust and dirt from getting in and minimizes outgassing or absorption of atmospheric gases. The complete cap assembly can be autoclaved.

Use low pressure [Connectors](#) (nuts) to secure Low Pressure [Tubing](#). Click [HERE](#) for more information about the plugs for this cap.

---

*Tubing is shown for display purposes only and is sold separately.*

---

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$48.15**

Printed On 2024/05/20