

## Catalog Number

59008-15P

**HPLC Column, MicroBee C8 Octyl**, 10um, 4.6mm ID x 150mm Length,  
125A. Cogent Brand. 1 EA.  
MOST POPULAR.



### Additional Info

This Cogent MicroBee C8™ Octyl Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is end capped along with traditional bonding of an 8 carbon straight chain ligand to irregular shaped silica. This makes this column an ideal Alternate in USP Methods that use the L7 Classification and Microbondapak™ or similar C8 columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is excellent for long chain fatty acids and will have less Retention than a C18 Column.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$775.46**

Printed On 2024/05/20