

## Catalog Number

59021-25P

**HPLC Column, Axis C18**, 5um, 4.6mm ID x 250mm Length, 100A.  
Cogent Brand. 1 EA.



### Additional Info

This Cogent Axis C18™ Stationary Phase has been a premium value and used by HPLC laboratories for many years and is a real legacy. It is double end capped (*fully*), with traditional bonding of an 18 carbon silane to the highest purity silica available, making this Column ideal as an Alternate in any USP Method that uses the L1 Classification.

With a 5 µm Particle, this column works well with higher flow rates and lower back pressures compared with smaller particles. Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids. Expect longer retention than with a C8 column.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$910.66**

Printed On 2024/05/20