

## Catalog Number

59300-02

**HPLC Column, MicroBee Phenyl, 10um, 3.9mm ID x 300mm Length, 125A. Cogent Brand. 1 EA.**  
**MOST POPULAR.**



### Additional Info

This Cogent MicroBee Phenyl™ Stationary Phase is part of the Similars™ product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A", 10um, 125A Stationary Phase is fully end capped and is manufactured with a 4 carbon ligand that has an attached Phenyl group. Bonded to irregular shaped silica, this makes it an ideal Alternate in USP Methods that use the L11 Classification and Microbondapak™ or similar Phenyl Columns.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$764.59**

Printed On 2024/05/19