

Catalog Number

59300-04

HPLC Column, MicroBee Amino, 10um, 3.9mm ID x 300mm Length,
125A. Cogent Brand. 1 EA.
MOST POPULAR.



Additional Info

This Cogent MicroBee Amino™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is not end capped and manufactured with traditional bonding of a 3 carbon chain ligand with a terminal Amino group, to irregular shaped silica. This Column is an ideal Alternate in USP Methods that use the L8 Classification and Microbondapak™ or similar Amino columns.

This Column is packed in and shipped in Normal Phase Solvents. This is one of our Most Popular MicroBee Amino Column.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$764.59

Printed On 2024/05/20