

Catalog Number

59350-01-RP

HPLC Column, MicroBee Amino, 10um, 3.9mm ID x 300mm Length, 125A. For Reversed Phase use. Cogent Brand. 1 EA.



Additional Info

This Cogent MicroBee Amino™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-A" Stationary Phase is not end capped and manufactured with traditional bonding of a 3 carbon chain ligand with a terminal Amino group, to irregular shaped silica. This Column is an ideal Alternate in USP Methods that use the L8 Classification and Microbondapak™ or similar Amino columns.

This Column is packed in and shipped in Reversed Phase Solvents.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$958.58

Printed On 2024/05/19