

## Catalog Number

62308-25P-4

**HPLC Column, Hyper Similar C8 Octyl**, 3 $\mu$ m, 4.0mm ID x 250mm Length, 150A. Cogent Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent Hyper Similar C8™ Octyl Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 8 carbon straight chain ligand to spherical shaped silica with very similar surface area to Hypersil® or Spherisorb® columns. This makes it an ideal Alternate in USP Methods that use the L7 Classification and Thermo Fisher Hypersil™ Columns or Similar C8 Columns.

The 3 $\mu$ m Particle in this column makes it more efficient than a similar 5 $\mu$ m column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure and methods that use longer columns and viscous solvents may exceed the pressure limits of your instrument. Used for hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is excellent for long chain fatty acids with less Retention than a C18 Column.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,036.93**

Printed On 2024/05/20