

## Catalog Number

62508-10P-4

**HPLC Column, Hyper Similar C8 Octyl**, 5um, 4.0mm ID x 100mm Length, 150A. Cogent Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent Hyper Similar C8™ Octyl Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 8 carbon straight chain ligand to spherical shaped silica with very similar surface area to Hypersil® or Spherisorb® columns. This makes it an ideal Alternate in USP Methods that use the L7 Classification and Thermo Fisher Hypersil™ Columns or Similar C8 Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octyl Phase is excellent for long chain fatty acids but expect less retention than with a C18 Column.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$785.53**

Printed On 2024/05/20