

Catalog Number

62518-15P-2

HPLC Column, Hyper Similars C18, 5um, 2.1mm ID x 150mm Length, 150A. Cogent Brand. 1 EA.



Additional Info

This Cogent Hyper Similars C18™ Stationary Phase is part of the Similars product line that has been a premium value used by HPLC laboratories for many years. The "Type-B" Stationary Phase is end capped along with traditional bonding of an 18 carbon straight chain ligand to spherical shaped silica with very similar surface area to Hypersil® or Spherisorb® columns. This makes it an ideal Alternate in USP Methods that use the L1 Classification and Thermo Fisher Hypersil™ Columns or Similar C18 Columns.

Used for hydrophobic or neutral compounds in Reversed Phase, this Octadecyl Phase is excellent for long chain fatty acids.

This narrow-bore column provides good run times but requires a pump that can precisely deliver, and a flow cell that is made for reduced flow rates with lower volumes. This column size is used for LCMS when more separation is desired

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70

Printed On 2024/05/20