

## Catalog Number

63008-10P-W

**HPLC Column, C8 Octyl, 5um, 4.6mm ID x 100mm Length, 100A with WCS hardware. Partisil Brand. 1 EA.**



### Additional Info

This Partisil C8™ Octyl, Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is modified with traditional bonding of a C8 ligand. This Column is written in many USP Methods; Partisil Silica used in older methods with the L7 Classification.

*Often used for hydrophobic or non polar compounds in Reversed Phase HPLC, this column is supplied with the Whatman Compression Sealing System RAC-II (WCS).*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,413.35**

Printed On 2024/05/20