

## Catalog Number

63087-25P-W

**HPLC Column, PAC, 5um, 4.6mm ID x 250mm Length, 100A with WCS hardware. Partisil Brand. 1 EA.**



### Additional Info

This Partisil PAC™ Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is irregular shaped, has a large surface area with high loading capacity and is not end capped and is modified with a ligand that contains both Amino and Cyano functional groups. This Column is written in many Methods that call for a Partisil Amino-Cyano (PAC) Column.

*This column is supplied with the Whatman Compression Sealing System RAC-II (WCS).*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,657.96**

Printed On 2024/05/20