

Catalog Number

65318-15P-2

HPLC Column, Zx Similar C18, 3um, 2.1mm ID x 150mm Length.
Cogent Brand. 1 EA.



Additional Info

This Cogent Zx Similar C18™ Stationary Phase is part of the Similar product line that has been a premium value used by HPLC laboratories for many years.

The "Type-B" Stationary Phase is end capped along with traditional bonding of an 18 carbon straight chain ligand to spherical shaped silica with very Similar surface area to Zorbax® columns. This makes it an ideal Alternate in USP Methods that use the L1 Classification and Agilent's Zorbax™ brand C18 HPLC Columns.

This narrow-bore column provides good run times but requires a pump that can precisely deliver, and a flow cell that is made for reduced flow rates with lower volumes. This column size is excellent, and a favorite for LCMS when additional resolution is desired.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,015.95

Printed On 2024/05/20