

Catalog Number

66000-C1

HPLC Column, Silica, 5um, 4.6mm ID x 125mm Length, 100A with WVS Cartridge hardware. Partisphere Brand. 1 EA.



Additional Info

This Partisphere Silica™ Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is spherically shaped and is not modified in any way. This Column is written in many USP Methods; Partisil Silica used in older methods with the L3 Classification.

Often used for hydrophilic or polar compounds in Normal Phase HPLC, this column is the cartridge style, Whatman Void Sealing System (WVS) and requires end fitting kit [66900-WEF](#) which is sold separately.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,032.84

Printed On 2024/05/20