

Catalog Number

66090-P25-W

HPLC Column, WAX Max-1, 5um, 4.6mm ID x 250mm Length, 100A with WVS Cartridge hardware. Partisphere Brand. 1 EA.



Additional Info

This Partisphere WAX™ Weak Anion Exchange Stationary Phase is the original, genuine Whatman / GE Healthcare product that has been a premium value used by HPLC laboratories for many years.

This "Type-A" Stationary Phase is spherically shaped, and is modified with a ligand that contains proprietary ionic groups. This Column is written in many USP Methods that call for a Partisphere WAX Column and the L23 Classification.

Often used for purification or analysis of biological compounds in Ion Exchange HPLC, this column is supplied with the cartridge style, Whatman Void Sealing System (WVS). This column requires the end fitting kit [66900-WEF](#) which is sold separately.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,291.03

Printed On 2024/05/20