

## Catalog Number

67091-05P-4

**HPLC Column, IEX WCX, 10um, 4.0mm ID x 50mm Length, 100A.**  
Cogent Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent Ion Exchange Column is a Weak Cation Exchanger (WCX has a negatively charged surface) Stationary Phase that has been a premium value and used by HPLC laboratories for many years. It is manufactured with ligands that include a carboxylic acid producing a functional charged range between pH 2.5 and 7.

We use only the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L9 Classification or for "base" analysis or purification.

With a 10µm Particle, this column works well with higher flow rates and lower back pressures for challenging IEX applications that include Cations (positively charged compounds).

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$691.26**

Printed On 2024/05/19