

## Catalog Number

68004-P21-PG1-3

**Prep Guard Column Cartridge, RP C4 Wide Pore, 10um, 300A.**  
21.2mm ID x 10mm Length. Cogent Brand. 1 EA.



### Additional Info

This Guard Column Cartridge is packed with 300Å, 10µm Cogent RP C4™ Stationary Phase. It is double end capped (*fully*), with traditional bonding of a 4 carbon ligand (Butyl) to the highest purity silica available. It is designed to protect your RP C4 Wide Pore, Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

*Use this Guard Column Cartridge with a 21.2mm ID Prep Column Holder. Click [HERE](#) for more information.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$447.95**

Printed On 2024/05/20