

Catalog Number

68008-SP-PG1

Semi Prep Guard Column Cartridge, RP C8 Octyl, 10um, 100A. 10mm ID x 10mm Length. Cogent Brand. 5/PK.



Additional Info

These Guard Column Cartridges are packed with 10µm, Cogent RP C8™ Octyl Stationary Phase. They are double end capped (*fully*), with traditional bonding of an 8 carbon silane to the highest purity silica available. They are designed to protect your RP C8 Semi Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, Stainless Steel Frits.

Use these Guard Column Cartridges with a 10mm ID Prep Column Holder. Click [HERE](#) for more information.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$895.91

Printed On 2024/05/19