

Catalog Number

68008-SP-PG1-3

Semi Prep Guard Column Cartridge, RP C8 Octyl, Wide Pore, 10um, 300A. 10mm ID x 10mm Length. Cogent Brand. 5/PK.



Additional Info

These Guard Column Cartridge are packed with 10µm, Cogent RP C8™ Octyl, Wide Pore, Stationary Phase. They are double end capped (*fully*), with traditional bonding of an 8 carbon silane to the highest purity silica available. They are designed to protect your RP C8 Wide Pore Semi Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, they have Stainless Steel Frits.

Use these Guard Column Cartridges with a 10mm ID Prep Column Holder. Click [HERE](#) for more information.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$895.91

Printed On 2024/05/19