

## Catalog Number

68022-P21-PG1

**Prep Guard Column Cartridge, RP PFP, 10um, 100A. 21.2mm ID x 10mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Guard Column Cartridge is packed with 10µm, Cogent RP PFP™ Pentafluorophenyl Stationary Phase. It is double end capped (*fully*), and bonded with a four carbon link to a benzene ring with 5 Fluorines on the highest purity 100A silica available.

It is designed to protect your RP PFP™ Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

*Use this Guard Column Cartridge with a 21.2mm ID Prep Column Holder. Click [HERE](#) for more information.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$447.95**

Printed On 2024/05/20