

Catalog Number

68024-P30-PG1

Prep Guard Column Cartridge, RP Phenyl, 10µm, 100A. 30mm ID x 10mm Length. Cogent Brand. 1 EA.



Additional Info

This Guard Column Cartridge is packed with 10µm, Cogent RP Phenyl™ Stationary Phase. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available. It is designed to protect your RP Phenyl Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

Use this Guard Column Cartridge with a 30mm ID Prep Column Holder click [HERE](#) for more information about the holder.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$537.54

Printed On 2024/05/19