

## Catalog Number

68024-SP-PG1

**Semi Prep Guard Column Cartridge, RP Phenyl, 10um, 100A. 10mm ID x 10mm Length. Cogent Brand. 5/PK.**



### Additional Info

These Guard Column Cartridges are packed with 10µm, Cogent RP Phenyl™ Stationary Phase. They are double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available. They are designed to protect your RP Phenyl Semi Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, they have SS Frits.

*Use these Guard Column Cartridges with a 10mm ID Prep Column Holder. Click [HERE](#) for more information.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$895.91**

Printed On 2024/05/19