

Catalog Number

68304-15P-2

HPLC Column, RP C4 Butyl, 3 μ m, 2.1mm ID x 150mm Length, 100A.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP C4™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 4 (*butyl*) carbon ligand to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L26 Classification.

With a 3 μ m Particle, this column is more efficient than a 5 μ m column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure. Used for very hydrophobic or neutral compounds in Reversed Phase, this Butyl Phase is excellent for short analyses of small molecules. Expect Less retention than with a C8 column.

A popular size, this narrow-bore column provides excellent data and requires a pump that can precisely deliver, and a flow cell that is made for reduced flow rates with lower volumes. This column size is excellent, and a favorite for LCMS.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,015.95

Printed On 2024/05/20