

Catalog Number

68322-10P-2

HPLC Column, RP PFP, 3um, 2.1mm ID x 100mm Length, 100A.

Cogent Brand. 1 EA.

MOST POPULAR.



Additional Info

This Very Popular Cogent RP PFP™ or RP Pentafluorophenyl Stationary Phase has been a premium value used by HPLC laboratories for many years and routinely called "F5". It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines; the highest purity silica available is used making this column ideal as an Alternate in USP

With a 3µm Particle, this column is more efficient than a 5µm column of the same phase and you will experience sharper peaks. However, due to the smaller particle size, there is an increased back pressure. This Column is often chosen for samples with halogenated or polar functional groups.

Packed and shipped in Reversed Phase Solvents.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70

Printed On 2024/05/20