

## Catalog Number

68504-P30-PG1-3

**Prep Guard Column Cartridge, RP C4 Wide Pore, 5µm, 300Å. 30mm ID x 10mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Guard Column Cartridge is packed with 300Å, 5µm Cogent RP C4™ Stationary Phase. It is double end capped (*fully*), with traditional bonding of a 4 carbon ligand (Butyl) to the highest purity silica available. It is designed to protect your RP C4 Wide Pore, Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

*Use this Guard Column Cartridge with a 30mm ID Prep Column Holder. Click [HERE](#) for more information.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$537.54**

Printed On 2024/05/19