

## Catalog Number

68508-P30-PG1

**Prep Guard Column Cartridge, RP C8 Octyl, 5um, 100A. 30mm ID x 10mm Length. Cogent Brand. 1 EA.**



### Additional Info

This Guard Column Cartridge is packed with 5µm, Cogent RP C8™ Octyl Stationary Phase. It is double end capped (*fully*), with traditional bonding of an 8 carbon silane to the highest purity silica available. It is designed to protect your RP C8 Octyl Prep Column from particulate and contaminants that will adsorb to it. Packed in a Stainless Steel Tube, it has SS Frits.

*Use this Guard Column Cartridge with a 30mm ID Prep Column Holder. Click [HERE](#) for more information.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$537.54**

Printed On 2024/05/19