

## Catalog Number

68508-P50-50-3M

**Process HPLC Column, RP C8 Octyl, Wide Pore, 5um, 300A. 50mm ID x 50mm Length. Cogent Brand. 1 EA.**  
*CUSTOM ORDER.*



### Additional Info

This Cogent RP C8™ Wide Pore, Stationary Phase has been a premium value, and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of an 8 (*Octyl*) carbon silane to the highest purity 300Å silica available.

The 5 µm Particle works well with higher flow rates and produces lower back pressures than smaller particles. The extra large diameter of this column provides high loading capacity which is great for purification of large compounds making it ideal for Process HPLC of many proteins & peptides and other compounds..

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$10,914.16**

Printed On 2024/05/19