

## Catalog Number

68522-P30-100

**Prep HPLC Column, RP PFP, 5um, 100A. 30mm ID x 100mm Length.**  
Cogent Brand. 1 EA.



### Additional Info

This Cogent RP PFP™ or RP Pentafluorophenyl Stationary Phase has been a premium value used by HPLC laboratories for many years and routinely called "F5". It is double end capped (*fully*) and bonded with a four carbon link to a benzene ring with 5 Fluorines on the highest purity 100A silica available.

With a 5µm Particle, this column works well with higher flow rates in preparative HPLC making this column ideal for purification of halogenated or polar samples.

*A specialized HPLC System is required for this column which is packed and shipped in Reversed Phase Solvents.*

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$7,369.23**

Printed On 2024/05/20