

Catalog Number

68524-10P

HPLC Column, RP Phenyl, 5um, 4.6mm ID x 100mm Length, 100A.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP Phenyl™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification or for new methods.

This size column is often used when shorter run times are needed but 50mm long columns do not provide enough resolution.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$785.53

Printed On 2024/05/19