

Catalog Number

68524-10P-2-3M

HPLC Column, RP Phenyl, Wide Pore, 5um, 2.1mm ID x 100mm
Length, 300A. Cogent Brand. 1 EA.



Additional Info

This Cogent RP Phenyl™ Wide Pore HPLC Column has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon chain with a terminal Phenyl Ring to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification for large molecules or in new method development.

Used to separate large Isomers (*positional*), Natural Products and Related Compounds and some Peptides in Reversed Phase HPLC.

This narrow-bore column provides good run times but requires a pump that can precisely deliver, and a flow cell that is made for reduced flow rates with lower volumes.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$848.36

Printed On 2024/05/19