

Catalog Number

68524-25P

HPLC Column, RP Phenyl, 5um, 4.6mm ID x 250mm Length, 100A.
Cogent Brand. 1 EA.



Additional Info

This Cogent RP Phenyl™ Stationary Phase has been a premium value and used by HPLC laboratories for many years. It is double end capped (*fully*), with traditional bonding of a 3 carbon silane with a terminal Phenyl Ring to the highest purity silica available, making this column ideal as an Alternate in USP Methods that use the L11 Classification or for new methods.

Used to separate Isomers (positional), Tocopherols, Flavonoids, Polynuclear Aromatics, Nitro-Aromatic Compounds, Active Pharmaceutical Ingredients (APIs), Natural Products and Related Compounds in Reversed Phase HPLC, these Phenyl Columns are an excellent choice.

This is Column is often used in USP Monographs.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$921.70

Printed On 2024/05/19