

Catalog Number

68540-P21-1003M

Prep HPLC Column, NP Amino, Wide Pore, 5um, 300A. 21.2mm ID x 100mm Length. Cogent Brand. 1 EA.



Additional Info

This Cogent NP Amino™ Wide Pore Stationary Phase has been a value and used by HPLC laboratories for many years. Not end capped, it is manufactured using a 3 carbon Ligand (propyl) with a terminal amino group, to the highest purity silica available. This makes this column ideal for preparative HPLC.

With a 5 µm Particle and 300A Pore Size, this column performs well for polar molecules with a negative charge such as complex Sugars and larger compounds with Carboxylic Acid Groups using Normal Phase or *HILIC* modes. These columns are packed and shipped in Normal Phase Solvents.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$3,684.61

Printed On 2024/05/20