

Catalog Number

69069-PG6

Guard Column cartridge, UDC Cholesterol, replacement, 4.0 mm ID x 10 mm long cartridge packed with 4 μ m, 100 Å phase. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your column. Packed in a Stainless Steel Tube, it has SS Frits. Labelled, and shipped in an individual, black, delrin case to protect and keep it from drying out. Use with [guard column holder](#) sold separately.

Sometimes, the individual cartridge is preferred over buying a pack of 5 when different "Guard Column" phases are preferred and pairing them with different analytical phases for unique separations.

Use this guard column cartridge with a 3.0, 4.0, and 4.6 mm ID analytical columns.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$244.79

Printed On 2024/05/19