

## Catalog Number

69220-10P-3

**HPLC Column, Phenyl Hydride 2.0**, 2.2 $\mu$ m, 3.0mm ID x 100mm Length, 120A, High Pressure. Cogent TYPE-C Brand. 1 EA. CUSTOM ORDER.



### Additional Info

This Cogent Phenyl Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for purification of compounds with aromatic rings that are difficult to retain or separate in Reversed Phase, including aromatic amines, small peptides and azides.

*This size column provides similar run times as a 4.6mm ID column, but uses much less mobile phase, lower flow rates and with less back pressure. No special equipment is needed to use this column.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,392.70**

Printed On 2024/05/20