

Catalog Number

69269-HG3

Guard Column Cartridge, UDC Cholesterol 2.o, Replacement
Cartridge, 2.0mm ID x 10mm Length, 2.2um, 120A. Hichrom style, in an individual black case. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your UDC Cholesterol 2.o™ Column. Packed in a Stainless Steel Tube, it has SS Frits. Labelled, and shipped in an individual, black, Delrin case to protect and keep it from drying out. Use with [Hichrom Style Guard Column Holder](#).

Sometimes, the individual cartridge is preferred over buying a pack of 5 when different "Guard Column" phases are preferred and pairing them with different analytical phases for unique separations.

Use this guard column cartridge with a 2.1mm ID analytical column packed with Undecanoic Cholesterol 2.o™ (UDC).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$275.19

Printed On 2024/05/19