

Catalog Number

70000-25P

HPLC Column, Diamond Hydride, 4um, 4.6mm ID x 250mm Length,
100A. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Diamond Hydride HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of polar compounds that are difficult to retain in Reversed Phase, including polar and non-polar peptides, amino acids and other biological compounds such as metabolites.

Chosen when more separation between closely related compounds or longer retention times is needed .

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,277.73

Printed On 2024/05/19